



QuEChERS

Quick, Easy, Cheap, Effective, Rugged
and Safe Method for Determining
Pesticide Residues



1



Fruits and Vegetables




Milk and Honey



3

Meat



4

Tobacco



5

QuEChERS

- Multiclass, multiresidue method (MRM)
- Primarily for the analysis of pesticides from hydrated (80-95% water content) matrices
- Multi-residue pesticide analysis of food and environmental samples can be problematic
- The complex sample matrix may contain abundant quantities of chlorophyll, lipids, sterols and other components that can interfere with good sample analysis
- Use of the QuEChERS method reduces these problems



The QuEChERS Method

- Consists of a liquid-liquid micro extraction
- Partitioning is aided by MgSO_4
- The preferred solvent is acetonitrile
- Acetonitrile provides extraction of the broadest range of organic compounds without co-extraction of large amounts of lipophilic material and is highly compatible with GC/MS and LC/MS applications showing the fewest interferences
- Followed by dispersive SPE sample clean-up to remove unwanted matrix materials

20 MINUTES!



Modifications to the Original Method

- Some modifications to the original QuEChERS method have been introduced to ensure efficient extraction of pH dependent compounds and to minimize degradation of base and acid labile pesticides
- Buffering with citrate salts has been introduced in the micro extraction to adjust the pH to a compromise value of 5 to 5.5, where most acid and base labile pesticides are sufficiently stabilized. To improve stability of base-labile compounds in the sample extracts, a small amount of formic acid is added to the final extract after cleanup
- Acidic pesticides are directly analyzed from the raw extract before PSA cleanup since they would be adsorbed and not released by the sorbent
- In another modification introduced by Schenck, a graphitized carbon black (GCB) PSA cartridge is used to remove plant pigments without the loss of planar compounds



Currently there are four variations of the QuEChERS method being used in western countries

- The original QuEChERS method. Introduced in 2003, uses sodium chloride to enhance extraction
- AOAC 2007.01. Uses sodium acetate as a buffer replacing sodium chloride and 1% acetic acid in acetonitrile
- The dual phase variation: This method variation introduces the use of a PSA and GCB cartridge to remove high levels of chlorophyll and plant sterols in the final extract without the loss of planar pesticides (polar aromatics) using an acetone:toluene solvent blend (3:1).
- The European version is similar to the AOAC method, except the extraction uses sodium chloride, sodium citrate dihydrate and disodium citrate sesquihydrate instead of sodium acetate.



Original QuEChERS

1. To **UCT** product **ECMSSC50CT** containing 6 grams anhydrous magnesium sulfate and 1.0 gram of sodium chloride in a 50 mL polypropylene centrifuge tube, add 15 ml of acetonitrile
2. Shake to mix contents
3. Add surrogate or internal standards if desired
4. Add 15 grams homogenized hydrated sample to the centrifuge tube
5. Shake for 1 minute
6. Centrifuge for 1 minute at 3700 rpm
7. Add a aliquot of the supernatant to the appropriate **UCT** dispersive solid-phase cleanup tube: (**CUMPSCB2CT**, **CUMPS2CT**, **CUMPSC18CT**, **ECMPSCB15CT**, or **ECMPSC1815CT**)
8. Shake for 1 minute
9. Centrifuge for 1 minute at 3700 rpm
10. Analyze extract





AOAC 2007.01

1. To **UCT** product **ECMSSA50CT** containing 6 grams anhydrous magnesium sulfate and 1.5 grams of anhydrous sodium acetate in a 50 mL polypropylene centrifuge tube, add 15 ml of 1% acetic acid in acetonitrile
2. Shake to mix contents
3. Add surrogate or internal standards if desired
4. Add 15 grams homogenized hydrated sample to the centrifuge tube
5. Shake for 1 minute
6. Centrifuge for 1 minute
7. Add an aliquot of the supernatant to the appropriate **UCT** dispersive clean-up tube: (**CUMPSCB2CT**, **CUMPS2CT**, **CUMPSC18CT**, **ECMPSCB15CT**, or **ECMPSC1815CT**)
8. Shake for 1 minute
9. Centrifuge for 1 minute at 3700 rpm
10. Analyze extract

11



The Dual Phase Variation

- Matrix plant pigments often interfere with analysis
- To reduce these interferences, graphitized carbon can be added to the dispersive solid-phase clean-up tubes
- Carbon however may result in a loss of planar (polar aromatic) pesticides
- Cleanup of plant pigments without loss of planar pesticides can be accomplished by using the **UCT** dual-phase cartridge containing GCB and PSA

12



Dual Phase Cartridge Clean-Up Procedure

1. Pre-rinse cartridge with 5 mL of toluene
2. Add an aliquot of the supernatant to the cartridge
3. Start collection
4. Elute with 6-12 mL of 3:1 acetone:toluene
5. Concentrate for GC/MS analysis or
6. Concentrate to dryness and reconstitute in mobile phase for LC analysis
7. Cartridge product selection used for this analysis:
ECPSACB6, **ECPSACB256** or ECPSACB506 depending upon sorbent mass required



Available Options from UCT

PART NUMBER	DESCRIPTION
CUMPSCB2CT	2mL micro-centrifuge tubes with 150mg Anhydrous Magnesium Sulfate, 50mg PSA & 50mg Carbon
CUMPS2CT	2mL micro-centrifuge tubes with 150mg Anhydrous Magnesium Sulfate, 50mg PSA
CUMPS18CT	2mL micro-centrifuge tubes with 150mg Anhydrous Magnesium Sulfate, 50mg PSA & 50mg endcapped C18
ECMPSCB15CT	15mL centrifuge tubes with 900mg Anhydrous Magnesium Sulfate, 300mg PSA & 150mg Carbon
ECMPSC1815CT	15mL centrifuge tubes with 900mg Anhydrous Magnesium Sulfate, 300mg PSA & 150mg endcapped C18
ECPSACB6	6mL columns with 400mg PSA on bottom, 200mg Graphitized Carbon on top, separated by a Teflon frit*
ECPSACB256	6mL columns with 250mg Graphitized Carbon on top, 500mg PSA on the bottom, separated with a Teflon frit*
ECPSACB506	6mL columns with 500mg Graphitized Carbon on top, 500mg PSA on the bottom, separated with a Teflon frit*
ECMSSA50CT	50mL PP centrifuge tube with 6g Anhydrous Magnesium Sulfate, 1.5g Anhydrous Sodium Acetate
ECMSSC50CT	50mL PP centrifuge tube with 4g Anhydrous Magnesium Sulfate, 1g Sodium Chloride
EUMIV50CT	50 mL PP centrifuge tube with 6g of magnesium sulfate anhydrous, 1.5g of sodium chloride, 1.5g of sodium citrate dihydrate, And 750 mg of disodium citrate sesquihydrate



Dispersive Products

Part Number	Contents and Use Description
CUMPSCB2CT	2mL micro-centrifuge tubes with 150mg Anhydrous Magnesium Sulfate, 50mg PSA & 50mg Carbon <i>A dispersive SPE product for removing pigments, polar organic acids, sugars and lipids. Will cause some loss of planar pesticides. Designed for use with a 2 mL aliquot of supernatant.</i>
CUMPS2CT	2mL micro-centrifuge tubes with 150mg Anhydrous Magnesium Sulfate, 50mg PSA <i>A dispersive SPE product for removing polar organic acids, some sugars and lipids. Designed for use with a 2 mL aliquot of supernatant.</i>
CUMPS18CT	2mL micro-centrifuge tubes with 150mg Anhydrous Magnesium Sulfate, 50mg PSA & 50mg endcapped C18 <i>A dispersive SPE product for removing polar organic acids, sterols, sugars and lipids. Designed for use with a 2 mL aliquot of supernatant.</i>

15



Dispersive Products Continued

Part Number	Contents and Use Description
ECMPSCB15CT	15mL centrifuge tubes with 900mg Anhydrous Magnesium Sulfate, 300mg PSA & 150mg carbon <i>A dispersive SPE product for removing pigments, polar organic acids, sugars and lipids. This product will cause the loss of planar pesticides. Designed for use with a 10 mL aliquot of supernatant.</i>
ECMPSC1815CT	15mL centrifuge tubes with 900mg Anhydrous Magnesium Sulfate, 300mg PSA & 150mg endcapped C18 <i>A dispersive SPE product for removing polar organic acids, sterols, some sugars and lipids from a 10 mL aliquot.</i>

16



Cartridge Products

Part Number	Contents and Use Description
ECPSACB6	6mL columns with 200mg Graphitized Carbon on top, 400mg PSA on bottom, separated by a Teflon frit <i>Used in the dual phase variation of QuEChERS, this product removes pigments, polar organic acids, some sugars and lipids from an aliquot of extract.</i>
ECPSACB256	6mL columns with 250mg Graphitized Carbon on top, 500mg PSA on the bottom, separated with a Teflon frit <i>Used for the same application as ECPSACB6 but with a different quantity of sorbents.</i>
ECPSACB506	6mL columns with 500mg Graphitized Carbon on top, 500mg PSA on the bottom, separated with a Teflon frit* <i>Used for the same application as ECPSACB6 but with a different quantity of sorbents.</i>

17



Extraction Products

ECMSSA50CT	50mL PP centrifuge tube with 6g Anhydrous Magnesium Sulfate, 1.5g Anhydrous Sodium Acetate <i>This product is listed in the AOAC for the extraction of pesticide residues using the QuEChERS method. It's primarily designed to preserve base sensitive compounds such as fungicides during extraction.</i>
ECMSSC50CT	50mL PP centrifuge tube with 4g anhydrous Magnesium Sulfate, 1.0 g Sodium Chloride <i>This product is used for the extraction of pesticide residues when base sensitive compounds are not an issue. The addition of sodium chloride improves extraction efficiency. This product also results in a cleaner extract.</i>
EUMIV50CT	50 mL polypropylene centrifuge tube with 6 g of anhydrous magnesium sulfate, 1.5 g of sodium chloride and 750 mg of disodium citrate sesquihydrate. <i>This product is the European version of ECMSSA50CT. It contains both a buffer and sodium chloride.</i>

18

The QuEChERS Method for Pesticide Residues




1) Shake sample with solvent and salts




2) Centrifuge for 1 min



3) Mix a portion with a sorbent



4) Centrifuge for 1 min



5) Analyze Pesticides


Add hydrated vegetation to salts and solvent SHAKE



Centrifuge



21



← MeCN Extract

← Plant Matrix


← Aqueous

Slide courtesy of Frank Schenck, FDA

22


Recoveries in the QuEChERS method

229 pesticides analyzed by GC-MS and LC-MS/MS



23

Recoveries of 15 Pesticides in Different Matrices




No differences were found vs. matrix for individual pesticides or concentration

24


Repeatability in the QuEChERS method

229 pesticides analyzed by GC-MS and LC-MS/MS




25

Why PSA CLEANUP (GC/FPD - Asparagus Extract)





Matrix components removed by SPE sorbents



Saito et al (2004)



Pesticide Analytes

GC amenable pesticides are capitalized. LC/MS/MS pesticides are not capitalized.
Analytes that can be analyzed under either instrument are underlined.

<u>acephate</u> *	acetamiprid	Acrinathrin	aldicarb	aldicarb sulfone
aldicarb sulfoxide	Aldrin	azaconazole	azamethiphos	<u>azinphos-methyl</u>
<u>azoxystrobin</u>	Bifenthrin	<u>bitertanol</u>	Bromopropylate	<u>bromuconazole</u>
Bupirimate	<u>buprofezin</u>	butocarboxim	butocarboxim sulfone	butocarboxim sulfoxide
Cadusafos	<u>carbaryl</u>	carbendazim	<u>carbofuran</u>	3-hydroxy-carbofuran
chlorbromuron	(α -, γ -)Chlordane	(α -, β -)Chlorfenvinphos	Chlorpropham	Chlorpyrifos
Chlorpyrifos-methyl	Chlorthalodimethyl	Chlorothalonil*	Chlorzolinate	clofentezine
Coumaphos	cycloxydim*	(λ -)Cyhalothrin	cymoxanil	Cypermethrin
<u>cyproconazole</u>	<u>cyprodinil</u>	(2,4'-4,4'-)DDE	(2,4'-4,4'-)DDT	Deltamethrin
demeton	demeton-O-sulfoxide	demeton-S-methyl	demeton-S-methyl sulfone	desmedipham
Diazinon	<u>dichlofluanid</u> *	Dichlorobenzophenone	<u>dichlorvos</u>	diclobutrazole
Dicloran	dicrotophos	Dieldrin	<u>Diethofencarb</u>	<u>difenoconazole</u>
Diflufenican	<u>dimethoate</u>	dimethomorph	<u>diniconazole</u>	Diphenyl
Diphenylamine	<u>disulfoton</u>	<u>disulfoton sulfone</u>	diuron	<u>dmsa</u>
<u>dmst</u>	dodemorph	α -Endosulfan	β -Endosulfan	Endosulfan sulfate
EPN	epoxiconazole	Esfenvalerate	etaconazole	ethiofencarb



ethiofencarb sulfoxide	Ethion	ethirimol	<u>Ethoprophos</u>	etofenprox
Etridiazole	Famoxadone	fenamiphos	fenamiphos sulfone	Fenarimol
Fenazaquin	fenbuconazole	fenhexamid*	Fenithrothion	fenoxycarb
Fenpiclonil	Fenpropathrin	Fenpropidine	fenpropimorph	fenpyroximate
<u>Fenthion</u>	<u>fenthion sulfoxide</u>	Fenvalerate	florasulam*	Flucythrinate I & II
Fludioxonil	flufenacet	Flufenconazole	flusilazole	Flutolanil
Fluvalinate	Fonophos	fosthiazate	Furalaxyl	furathiocarb
furmecyclox	Heptachlor	Heptachlor epoxide	Heptenophos	Hexachlorobenzene
hexaconazole	hexythiazox	imazalil	imidacloprid	Iprodione
iprovalicarb	isoprothiolane	isoxathion	kresoxim-methyl	Lindane
linuron	<u>Malathion</u>	<u>malathion oxon</u>	Mecarbam	mephosfolan
Mepronil	Metaxyl	metconazole	methamidophos*	Methidathion
<u>methiocarb</u>	methiocarb sulfone*	methiocarb sulfoxide	methomyl	methomyl-oxime
metobromuron	metoxuron	Mepanipyrim	Mevinphos	<u>monocrotophos</u>
monolinuron	<u>myclobutanil</u>	nuarimol	Ofurace	<u>omethoate</u>
oxadixyl	oxamyl	oxamyl-oxime	oxydemeton-methyl	paclobutrazole



Parathion	Parathion-methyl	<u>penconazole</u>	<u>pencycuron</u>	cis- Permethrin
trans-Permethrin	phenmedipham	o-Phenylphenol	Phorate	<u>phorate sulfone</u>
Phosalone	Phosmet	Phosmet-oxon	phosphamidon	Phthalimide
<u>picoxystrobin</u>	Piperonyl butoxide	<u>pirimicarb</u>	<u>pirimicarb-desmethyl</u>	Pirimiphos-methyl
prochloraz	Procymidone	<u>profenofos</u>	Prometryn	Propargite
Propham	<u>propiconazole</u>	<u>propoxur</u>	Propyzamide	Prothiofos
pymetrozine*	Pyrazophos	<u>pyridaben</u>	<u>pyridaphenthion</u>	<u>pyrifenoxy</u>
pyrimethanil	Pyriproxyfen	Quinalphos	Quinoxifen	Quintozene
sethoxydim*	spinosad	<u>spiroxamine</u>	<u>tebuconazole</u>	tebufenozide
<u>Tebufenpyrad</u>	<u>tetraconazole</u>	Tetradifon	Tetrahydrophthalimide	Terbufos
Terbufos sulfone	thiabendazole	thiacloprid	thiamethoxam	thiodicarb
thiofanox	thiofanox sulfone	thiofanox sulfoxide	thiometon	thiometon sulfone
thiometon sulfoxide	thiophanate-methyl	Tolclofos-methyl	<u>tolylfluanid*</u>	<u>triadimefon</u>
<u>triadimenol</u>	Triazophos	trichlorfon	tricyclazole	tridemorph
<u>trifloxystrobin</u>	<u>trifluminazole</u>	Trifluralin	<u>Triphenylphosphate</u>	vamidothion
vamidothion sulfone	vamidothion sulfoxide	Vinclozolin		




Why buy it from UCT?

Since many laboratories assemble their own clean-up products for this analysis, a comparison study was conducted at the USDA ARS Eastern Regional Research Center to determine if commercially prepared QuEChERS products could be as effective as products prepared in the users lab. Bulk magnesium sulfate anhydrous, primary secondary amine (PSA) and endcapped C18 were purchased from a commercial source and compared to the commercial product assembled using the same lot of bulk sorbents. The ratio of magnesium sulfate, PSA and C18 was 3:1:1. The clean-up products were tested on extracts of milk, honey and soybean. Efficacy of clean-up was determined by GC/MS analysis and compared the number of peaks above threshold values. Results clearly showed that the commercially prepared product provided fewer interferences than the product prepared in the lab in all three matrices. The extra peaks observed in the lab prepared product were probably caused by contamination from the lab environment. The commercially assembled product was prepared under controlled conditions minimizing potential contamination. These results, coupled with the obvious time and labor savings for assembly, indicate that commercially available QuEChERS products are preferable to products made "in-house".


31



Honey Extract Cleaned with "In House" Product



Honey Extract Cleaned with Commercial Product



Chromatograms for soybean and milk products showed similar improved clean-up when using commercial vs. "in-house" prepared products.


GC Peak Count

Above Threshold Value

Matrix	HONEY		SOYBEAN		MILK	
	UCT	IN-HOUSE	UCT	IN-HOUSE	UCT	IN-HOUSE
Replicate	# of peaks	# of peaks	# of peaks	# of peaks	# of peaks	# of peaks
1	7	20	7	17	43	91
2	9	12	8	15	49	103
3	7	21	5	20	52	108
4	8	24	2	12	43	121
5	5	18	6	8	46	117
6	5	22	2	13	45	104
7	8	8	7	11	49	89
8	4	13	4	10	103	117
9	5	18	4	7	107	127
10	8	12	3	9	106	127
11	6	15	2	31	116	120
12	6	12	8	28	126	118
13	6	19	6	35	104	119
14	6	21	4	51	106	108
15	5	20	4	43	100	118
16	4	14	7	43	109	113
Average	6	17	5	22	81	113

33

Summary Comparison of Sample Cleanliness
based upon GC peak count



34

Conclusions

- GC peak counts were significantly lower with the commercial QuEChERS product, indicating cleaner extracts
- Variability of peak counts on replicate samples were significantly lower using commercial QuEChERS products
- Savings in preparation costs and reduced product variability when using commercial QuEChERS products

References

M. Anastassiades, S.J. Lehotay, D. Stainbahr, and F. Schenck, (2003) J. AOAC Internat., 86, 412-431.

Schenck, F.J. and Vega, V. (2001) The Elution/Retention Characteristics of Pesticides on Graphitized Carbon and Amino Solid Phase Extraction Columns. LIB# 4244

S.J. Lehotay, "Quick, Easy, Cheap, Effective, Rugged and Safe (QuEChERS) Approach for Determining Pesticide Residues" in Methods in Biotechnology, Vol. 19, Pesticide Protocols, Martinez Vidal, J.L. and Garrido Frenich, A. (Eds.) Humana Press, Totowa, NJ, 2005, pp. 239-261.

F. J. Schenck, Victor Vega, "The Elution/Retention Characteristics of Pesticides on Graphitized Carbon and Amino Solid Phase Extraction Columns" Southeastern Regional Laboratory Pesticides; April 2001 Laboratory Information Bulletin

References

S.J. Lehotay, A.de Kok, M. Hiemstra, and P. Van Bodegraven, "Validation of a Fast and Easy Method for the Determination of 229 Pesticide Residues in Fruits and Vegetables Using Gas and Liquid Chromatography and Mass Spectrometric Detection", J. AOAC Int., 88, 595-614 (2005).

S.J. Lehotay, K. Mastovska, and A.R. Lightfield, "Use of buffering to improve results of problematic pesticides in a fast and easy method for residue analysis of fruits and vegetables" J. AOAC Int., 88, 615-629 (2005).

S.J. Lehotay, K. Mastovska, and S.J. Yun, "Evaluation of Two Fast and Easy Methods for Pesticide Residue Analysis in Fatty Food Matrices" J. AOAC Int., 88, 630-638 (2005)

F.J. Schenck and J.E. Hobbs, "Evaluation of the Quick, Easy, Cheap, Effective, Rugged and Safe (QuEChERS) approach to pesticide residue analysis" Bull. Environ. Contam. Toxicol., 73, 24-30 (2004)

Current Trends in QuEChERS Research

Multiresidue analytical method using dispersive solid-phase extraction and gas chromatography/ion trap mass spectrometry to determine pharmaceuticals in whole blood

Florian Plossl, Martin Giera, Franz Bracher
Department Pharmazie – Zentrum für Pharmaforschung, Ludwig-Maximilians-Universität München,
Butenandtstrasse 5-13, 81377 Munich, Germany

QuEChERS for Anthelmintics

New Method for the Analysis of Anthelmintics in Animal Tissues


**Brian Kinsella^{1,2,3}, Steven Lehotay³, Katerina Mastovska³,
Alan Lightfield³, Martin Danaher¹, Ambrose Furey²**

¹Teagasc

²Cork Institute of Technology

³Eastern Regional Research Center,
USDA Agricultural Research Service

brian.kinsella@ars.usda.gov





Extraction and Clean-up

1. Add 10 mL acetonitrile to **UCT** product **ECMSSC50CT**
2. Shake
3. Add 10 g. hydrated sample
4. Add internal standard (Cyprodinil + 2,4D)
5. Shake for 1 minute
6. Centrifuge for 5 minutes
7. Clean-up by adding 1 mL of supernatant to a 2 mL centrifuge containing 50 mg **UCT EEC18** sorbent and 150 mg **UCT MgSO₄**
8. Mix for 1 minute
9. Centrifuge for 1 minute



Analysis

- Place 0.5 mL into auto sampler vial
- Add QC spike (TPP)
- Inject onto LC-MS/MS (ESI+ and ESI-)

39 Flukicides/Anthelmintics

ESI +

Triphenyl Phosphate (IS)
Cyprodinil (QC spike)

Abamectin
Doramectin
Emamectin
Eprinomectin
Moxidectin
Ivermectin
Selamectin

Dichlorvos
Coumaphos
Coumaphos-Oxon
Haloxon
Morantel
Levamisole

Albendazole
Albendazole-Sulfoxide
Albendazole-Sulfone
Albendazole-Amino-Sulfone
C mebendazole
Flubendazole
Amino-Flubendazole
Hydroxy-Flubendazole
Fenbendazole
Fenbendazole-Sulfone
Oxfendazole
Mebendazole
Amino-Mebendazole
Hydroxy-Mebendazole
Oxibendazole
Thiabendazole
5-Hydroxy-Thiabendazole
Triclabendazole
Triclabendazole-Sulfoxide

ESI -

2,4-D (IS)

Bithionol
Clorsulon
Closantel
Niclosamide
Nitroxynil
Oxyclozanide
Rafoxanide
Triclabendazole

Slide adapted from Brian Kinsella, Teagasc

Removal of Co-extractives


Percentage of Co-extractives Removed

	PSA	C18	PSA+C18
Liver	81%	36%	89%
Milk	10%	73%	90%

- C₁₈ removes lipophilic compounds
- PSA removes organic acids
- PSA + C₁₈ removed highest % of co-extractives
- But PSA + C₁₈ also removed analytes
- C₁₈ alone provides sufficient clean-up and gives high recoveries


Slide adapted from Brian Kinsella, Teagasc

Recovery Experiments




Slide adapted from Brian Kinsella, Teagasc

Recovery Experiments



Slide adapted from Brian Kinsella, Teagasc

Recoveries of the 39 Analytes



Slide adapted from Brian Kinsella, Teagasc

Vet. Drug Analysis in Animal Tissue






Streamlined sample preparation method for analysis of many antibiotics in beef kidney/juice or serum

- Weigh out 1 g of homogenized beef kidney sample, or kidney juice or serum in a 50 mL FEP (fluorinated ethylene propylene) tube (a disposable polypropylene Corning tube can be used instead).
- Add 100 μ L of 1 μ g/mL composite internal standard solution of 13 C-sulfamethazine (serving for volume changes compensation), penicillin-V and cefadroxil (a penicillin and cephalosporin, respectively, serving for method performance control) in water.
- Add 2 mL water and 8 mL acetonitrile.
- Vortex briefly, shake for 5 min.
- Centrifuge at 3450 rcf for 5 min.
- Decant the supernatant into a 50 mL tube with 500 mg **UCT EEC18 bulk sorbent**.
- Vortex briefly, shake for 30 s.
- Centrifuge at 3450 rcf for 1 min.
- Place 5 mL aliquot of the supernatant into a graduated tube.
- Evaporate down to < 1 mL.
- Make up the volume to 1 mL with water.
- Transfer the extract into Whatman Uni-Prep vials (PVDF, 0.45 μ m) and compress the plunger to filter the extracts, which are then ready for LC-MS/MS analysis.



Recoveries of selected sulfonamide, fluoroquinolone, macrolide, and tetracycline antibiotics (50 ng/g, beef kidney)




Slide adapted from Kate Mastovska, USDA-ARS

Recoveries of β -Lactam Antibiotics

(50 ng/g, beef kidney)

■ normalized to Pen V

■ not normalized



DCCD = desfurylethiofur cysteine disulfide

Slide by Kate Mastovska, USDA-ARS


LC-MS/MS results for 121 veterinary drugs in kidney juice

spiked at 10, 25, 50, and 100 ng/mL

n = 16 (1 spike each day for 16 days)

80-120% avg recov (10-20% RSD):	35 analytes
70-120% avg recov (21-30% RSD):	52
50-69% avg recov (20-30% RSD):	9
55-130% avg recov (31-42% RSD):	17
20-35% avg recov (30-45% RSD):	6
Good results but degradation issue:	2
(cephapirin => desacetyl cephapirin)	

QuEChERS without Rugged is:



Curtain plate after injection of 25 samples with extracts from raisins without cleanup

Picture from Lutz Alder, BfR, Germany

53

Don't Suffer in Silence!
This man's job is to help you.



Contact him with any questions or problems.

dshelly@unitedchem.com

54